

PATENT APPLICATION

**PULSE OXIMETER WITH ALTERNATE HEART-RATE
DETERMINATION**

Inventor: **Clark R. Baker, Jr.**, a citizen of the United States, residing at
18493 Magee Way
Castro Valley, California 94546

Assignee: **Nellcor Puritan Bennett Incorporated**
4280 Hacienda Drive
Pleasanton, California 94588,
a corporation of Delaware

Entity: Large

**PULSE OXIMETER WITH ALTERNATE HEART-RATE
DETERMINATION**

CROSS-REFERENCES TO RELATED APPLICATIONS

5 [0001] NOT APPLICABLE

**STATEMENT AS TO RIGHTS TO INVENTIONS MADE UNDER
FEDERALLY SPONSORED RESEARCH OR DEVELOPMENT**

[0002] NOT APPLICABLE

10

**REFERENCE TO A "SEQUENCE LISTING," A TABLE, OR A COMPUTER
PROGRAM LISTING APPENDIX SUBMITTED ON A COMPACT DISK**

[0003] NOT APPLICABLE

15

BACKGROUND OF THE INVENTION

[0004] The present invention relates to oximeters, and in particular to determining a pulse rate by multiple mechanisms in a detected waveform from a pulse oximeter.

20

[0005] Pulse oximetry is typically used to measure various blood chemistry characteristics including, but not limited to, the blood-oxygen saturation of hemoglobin in arterial blood, the volume of individual blood pulsations supplying the tissue, and the rate of blood pulsations corresponding to each heartbeat of a patient. Measurement of these characteristics has been accomplished by use of a non-invasive sensor which scatters light through a portion of the patient's tissue where blood perfuses the tissue, and photoelectrically senses the absorption of light in such tissue. The amount of light absorbed at various wavelengths is then used to calculate the amount of blood constituent being measured.

25

30

[0006] The light scattered through the tissue is selected to be of one or more wavelengths that are absorbed by the blood in an amount representative of the amount of the blood constituent present in the blood. The amount of transmitted light scattered through the tissue will vary in accordance with the changing amount of blood constituent in the tissue and the related light absorption. For measuring blood oxygen level, such sensors have typically been

provided with a light source that is adapted to generate light of at least two different wavelengths, and with photodetectors sensitive to both of those wavelengths, in accordance with known techniques for measuring blood oxygen saturation.

[0007] Known non-invasive sensors include devices that are secured to a portion of the body, such as a finger, an ear or the scalp. In animals and humans, the tissue of these body portions is perfused with blood and the tissue surface is readily accessible to the sensor.

[0008] Patents No. 6,083,172, No. 5,853,364 and No. 6,411,833 show multiple methods of calculating a pulse rate in a pulse oximeter, with a "best rate" module which arbitrates between the pulse rate calculations to select a best rate based on confidence levels associated with each. The confidence levels are calculated using various metrics to determine the reliability of the different pulse rate calculations. Also, Patent No. 5,524,631 shows a fetal heart rate monitor that uses multiple parallel filter paths to identify the fetal heart rate, and uses a figure of merit operation to weight the different heart rate estimates.

[0009] N-100. The N-100 technology, dating to around 1985, accepted or rejected pulses based on pulse history of the size of pulses, pulse shape, expected time to occur (frequency) and ratio of R/IR.

[0010] In particular, the N-100 found pulses by looking for a signal maximum, followed by a point of maximum negative slope, then a minimum. The processing was done in a state machine referred to as "munch." Each maximum was not qualified until the signal passed below a noise threshold, referred to as a noise gate. This acted as an adaptive filter since the noise gate level was set by feedback from a subsequent processing step to adapt to different expected signal amplitudes. The pulses are then accepted or rejected in a "Level3" process which was a filter which adapts to changing signals by comparing the amplitude, period and ratio-of-ratios (ratio of Red to IR, with Red and IR being expressed as a ratio of AC to DC) of a new pulse to the mean of values in a history buffer, then determining if the difference is within a confidence level. If the new pulse was accepted, the history buffer was updated with the values for the new pulse. The level3 process acted as an adaptive bandpass filter with center-frequency and bandwidth (confidence limits) being adapted by feedback from the output of the filter.

[0011] N-200. The N-200 improved on the N-100 since it could be synchronized with an ECG, and included ECG filtering. The N-200 also added interpolation to compensate for baseline shift between the time of measuring the pulse maximum and minimum. The N-200

included other filtering features as well, such as a "boxcar" filter which computed the mean of a varying number of signal samples.

[0012] The N-200, after various filtering and scaling steps, applies the digitized signals to a "boxcar" filter, which computes the mean of N samples, where N is set by feedback from a subsequent processing step according to the filtered heart rate. New samples are averaged into the boxcar filter, while the oldest samples are dropped. The boxcar length (N) is used to set three parameters: a pulse threshold, absolute minimum pulse and small pulse. An ensemble-averaging (a.k.a "slider") filter then produces a weighted average of the new samples and the previous ensemble-averaged sample from one pulse-period earlier. The samples are then passed to a "munch" state machine and a noise gate, like the N-100. An interpolation feature is added to the N-100 process, to compensate for changes in the baseline level. Since the minimum and maximum occur at different times, a changing baseline may increase or decrease the minimum and not the maximum, or vice-versa.

[0013] "Ensemble averaging" is an integral part of C-Lock, which is NELLCOR's trademark for the process of averaging samples from multiple pulses together to form a composite pulse. This process is also known as "cardiac-gated averaging." It requires a "trigger" event to mark the start of each pulse.

BRIEF SUMMARY OF THE INVENTION

[0014] The present invention is a pulse oximeter which determines multiple heart rates, and selects between them based on the metrics of only one of the heart rate calculations. A primary heart rate calculation method is selected, and is used unless its metrics indicate questionable accuracy, in which case an alternative rate calculation is available and is used instead.

[0015] In one embodiment, the primary heart rate calculation method does not use an ensemble averaged waveform, while the alternative heart rate calculation does use an ensemble averaged waveform. The alternative heart rate calculation is used if the primary calculation has disqualified its most recently detected pulse.

BRIEF DESCRIPTION OF THE DRAWINGS

[0016] Fig. 1 is a block diagram of an oximetry system incorporating an embodiment of the invention.

[0017] Fig. 2 is a diagram of the software processing blocks of an oximeter including an embodiment of the present invention.

[0018] Fig. 3 is a context diagram of the pulse rate calculation subsystem.

5

DETAILED DESCRIPTION OF THE INVENTION

[0019] Fig. 1 illustrates an embodiment of an oximetry system incorporating the present invention. A sensor 10 includes red and infrared LEDs and a photodetector. These are connected by a cable 12 to a board 14. LED drive current is provided by an LED drive interface 16. The received photocurrent from the sensor is provided to an I-V interface 18.

10 The IR and red voltages are then provided to a sigma-delta interface 20 incorporating the present invention. The output of sigma-delta interface 20 is provided to a microcontroller 22 which includes a 10-bit A/D converter. Controller 22 includes flash memory for a program, and EEPROM memory for data. The processor also includes a controller chip 24 connected to a flash memory 26. Finally, a clock 28 is used and an interface 30 to a digital calibration
15 in the sensor 10 is provided. A separate host 32 receives the processed information, as well as receiving an analog signal on a line 34 for providing an analog display.

[0020] Design Summary The design of the present invention is intended to deal with unwanted noise. Signal metrics are measured and used to determine filter weighting. Signal metrics are things that indicate if a pulse is likely a plethysmograph or noise, such as
20 frequency (is it in the range of a human heart rate), shape (is it shaped like a heart pulse), rise time, etc. A similar technique was used in the Nellcor N200, described in the background of this application. The new design adds a number of different features and variations, such as the use of two ensemble averagers as claimed in the present invention.

[0021] Details of the architecture are shown in the diagram of Fig. 2. This design
25 calculates both the oxygen saturation, and the pulse rate, which are described separately below.

[0022] I. Oxygen Saturation Calculation.

[0023] A. Signal Conditioning - The digitized red and IR signals are received and are conditioned in this block by (1) taking the 1st derivative to get rid of baseline shift, (2) low
30 pass filtering with fixed coefficients, and (3) dividing by a DC value to preserve the ratio. The function of the Signal Conditioning subsystem is to emphasize the higher frequencies

that occur in the human plethysmograph and to attenuate low frequencies in which motion artifact is usually concentrated. The Signal Conditioning subsystem selects its filter coefficients (wide or narrow band) based on hardware characteristics identified during initialization.

5 [0024] Inputs - digitized red and IR signals

[0025] Outputs - Pre-processed red and IR signals

[0026] B. Pulse Identification and Qualification - The low pass filtered and digitized red and IR signals are provided to this block to identify pulses, and qualify them as likely arterial pulses. This is done using a pre-trained neural net, and is primarily done on the IR signal.

10 The pulse is identified by examining its amplitude, shape and frequency, just as was done in the Nellcor N-100. An input to this block is the average pulse period from block D. This function is similar to the N-100, which changed the upfront qualification using the pulse rate. The output indicates the degree of arrhythmia and individual pulse quality.

[0027] Inputs - (1) Pre-processed red and IR signals, (2) Ave. pulse period, (3)

15 Lowpass Waveforms from the low pass filter.

[0028] Outputs - (1) Degree of arrhythmia, (2) pulse amplitude variations, (3) individual pulse quality, (4) Pulse beep notification, (5) qualified pulse periods and age.

[0029] C. Compute Signal Quality Metrics - This block determines the pulse shape (derivative skew), period variability, pulse amplitude and variability, Ratio of Ratios
20 variability, and frequency content relative to pulse rate.

[0030] Inputs - (1) raw digitized red and IR signals, (2) degree of arrhythmia, individual pulse quality, pulse amplitude variation (3) pre-processed red and IR signals, (4) average pulse period.

[0031] Outputs - (1) Lowpass and ensemble averaging filter weights, (2) metrics for
25 sensor off detector, (3) Normalized Pre-processed waveforms, (4) percent modulation.

[0032] D. Average Pulse Periods. This block calculates the average pulse period from the pulses received.

[0033] Inputs - Qualified pulse periods and age.

[0034] Outputs - Average pulse period.

[0035] E1. Lowpass Filter and Ensemble Averaging - Block E1 low pass filters and ensemble averages the signal conditioned by block A, and normalized by block C, for the pulse rate identification. The weights for the low pass filter are determined by the Signal Metrics block C. The signal is also ensemble averaged (this attenuates frequencies other than those of interest near the pulse rate and its harmonics), with the ensemble averaging filter weights also determined by Signal Metrics block C. Less weight is assigned if the signal is flagged as degraded. More weight is assigned if the signal is flagged as arrhythmic because ensemble-averaging is not appropriate during arrhythmia. Red and IR are processed separately, but with the same filtering weights. The filtering is delayed approximately one second to allow the signal metrics to be calculated first.

[0036] The filters use continuously variable weights. If samples are not to be ensemble-averaged, then the weighting for the previous filtered samples is set to zero in the weighted average, and the new samples are still processed through the code. This block tracks the age of the signal - the accumulated amount of filtering (sum of response times and delays in processing). Too old a result will be flagged (if good pulses haven't been detected for awhile).

[0037] Inputs - (1) normalized pre-processed red and IR signals, (2) average pulse period, (3) low pass filter weights and ensemble averaging filter weights, (4) ECG triggers, if available, (5) IR fundamental, for zero-crossing triggers.

[0038] Outputs - (1) filtered red and IR signals, (2) age.

[0039] F. Estimate Filtered Waveform Correlation and Calculate Averaging Weight - this uses a noise metric similar to that used in the N100 and N200 described above, and doesn't use feedback. The variable weighting for the filter is controlled by the ratio-of-ratios variance. The effect of this variable-weight filtering is that the ratio-of-ratios changes slowly as artifact increases and changes quickly as artifact decreases. The subsystem has two response modes. Filtering in the Fast Mode targets an age metric of 3 seconds. The target age is 5 seconds in Normal Mode. In Fast Mode, the minimum weighting of the current value is clipped at a higher level. In other words, a low weight is assigned to the newest ratio-of-ratios calculation if there is noise present, and a high weight if no noise is present.

[0040] Inputs - (1) filtered red and IR signals and age, (2) calibration coefficients, (3) response mode (user speed settings).

[0041] Outputs - averaging weight for ratio-of-ratios calculation.

[0042] H. Calculate Saturation - Saturation is calculated using an algorithm with the calibration coefficients and averaged ratio of ratios.

[0043] Inputs - (1) Averaged Ratio-of-Ratios, (2) calibration coefficients.

5 [0044] Outputs - Saturation.

[0045] II. Pulse Rate Calculation.

[0046] E2. Lowpass Filter and Ensemble Averaging - Block E2 low pass filters and ensemble averages the signal conditioned by block A, for the pulse rate identification. The weights for the low pass filter are determined by the Signal Metrics block C. The signal is
10 also ensemble averaged (this attenuates frequencies other than those of interest near the pulse rate and its harmonics), with the ensemble averaging filter weights also determined by Signal Metrics block C. Less weight is assigned if the signal is flagged as degraded. More weight is assigned if the signal is flagged as arrhythmic since filtering is not appropriate during arrhythmia. Red and IR are processed separately. The process of this block is delayed
15 approximately one second to allow the signal metrics to be calculated first.

[0047] The filters use continuously variable weights. If samples are not to be ensemble-averaged, then the weighting for the previous filtered samples is set to zero in the weighted average, and the new samples are still processed through the code. This block tracks the age of the signal - the accumulated amount of filtering (sum of response times and delays in
20 processing). Too old a result will be flagged (if good pulses haven't been detected for awhile).

[0048] Inputs - (1) pre-processed red and IR signals, (2) average pulse period, (3) Lowpass filter weights and ensemble averaging filter weights, (4) ECG triggers, if available, (5) IR fundamental, for zero-crossing triggers.

25 [0049] Outputs - (1) filtered red and IR signals, (2) age.

[0050] I. Filtered Pulse Identification and Qualification - This block identifies and qualifies pulse periods from the filtered waveforms, and its results are used only when a pulse is disqualified by block B.

[0051] Inputs - (1) filtered red and IR signals and age, (2) average pulse period, (3)
30 hardware ID or noise floor, (4) kind of sensor.

[0052] Outputs - qualified pulse periods and age.

[0053] J. Average Pulse Periods and Calculate Pulse Rate - This block calculates the pulse rate and average pulse period.

[0054] Inputs - Qualified pulse periods and age

5 [0055] Outputs - (1) average pulse period, (2) pulse rate.

[0056] III. Venous Pulsation

[0057] K. Detect Venous Pulsation - Block K receives as inputs the pre-processed red and IR signal and age from Block A, and pulse rate and provides an indication of venous pulsation as an output. This subsystem produces an IR fundamental waveform in the time
10 domain using a single-tooth comb filter which is output to the Ensemble Averaging filters.

[0058] Inputs - (1) filtered red and IR signals and age, (2) pulse rate.

[0059] Outputs - Venous Pulsation Indication, IR fundamental

[0060] IV. Sensor Off

[0061] L. Detect Sensor-Off and Loss of Pulse Amplitude - The Pulse Lost and Sensor Off
15 Detection subsystem uses a pre-trained neural net to determine whether the sensor is off the patient. The inputs to the neural net are metrics that quantify several aspects of the behavior of the IR and Red values over the last several seconds. Samples are ignored by many of the oximetry algorithm's subsystems while the Signal State is not either Pulse Present or Sensor Maybe Off. The values of the Signal State variable are: "Pulse Present, Disconnect, Pulse
20 Lost, Sensor Maybe Off, and Sensor Off."

[0062] Inputs - (1) metrics, (2) front-end servo settings and ID

[0063] Outputs - Signal state including sensor-off indication

Pulse Rate Calculation subsystem

[0064] The subsystem averages qualified pulse periods from the Pulse Identification and
25 Qualification subsystem. It outputs the average period and the corresponding pulse rate.

[0065] The oximetry algorithm contains two instances of this subsystem. The first instance receives input from the Pulse Identification and Qualification instance whose input waveform have been processed by the Signal Conditioning subsystem, then lowpass filtered, but not

ensemble averaged, by the Ensemble Averaging subsystem. The second instance of the Pulse Rate Calculation subsystem receives input from two instances of the Pulse Identification and Qualification subsystem, the one described above and a second instance that receives input that has been ensemble averaged.

5 Selection of Pulse Period Source

[0066] One instance of the subsystem receives qualified pulse periods from two sources. The subsystem selects which of these two sources to use for its pulse rate calculation based solely on analysis of only one source, the “primary” source. The oximetry algorithm designates the Pulse Identification and Qualification instance that does NOT receive
10 ensemble-averaged waveforms as the primary source, and designates the other Pulse Identification and Qualification instance as the “alternate” source of qualified pulse periods. Qualified pulse periods from the alternate source are only used if the most recent pulse from the primary source was rejected. When a qualified pulse period is received from the primary
15 source, it is always used to update the pulse-rate calculation, and will prevent qualified pulse periods from the alternate source from being used until the primary source once again rejects a pulse period.

Calculation of Average Pulse Period and Pulse-Rate Estimate

[0067] When the subsystem uses a *Qualified_Pulse_Period*, it updates its average pulse period, *Avg_Period*, using a pulse-based, variable-weight IIR filter, then computes its *Rate*
20 output from *Avg_Period*. The steps for this filtering operation are:

1. $r_i = (60 / \Delta t) / \text{Qualified_Pulse_Period}$
2. $k = \text{Consecutive_Qualified} / \max(|r_i - r_{i-1}|, |r_{i-1} - r_{i-2}|, |r_{i-2} - r_{i-3}|, 1.0)$
3. $x = \text{bound}(\min(\text{Avg_Period}_{i-1}, \text{Qualified_Pulse_Period}), \frac{3}{4} \text{ seconds}, 2 \text{ seconds}) / 7 \text{ seconds}$
4. If $\text{Rate_Age} > 10 \text{ seconds}$, $x = \min(x * \text{Rate_Age} / 10 \text{ seconds}, 0.3)$
- 25 5. $k = \max(1 / \text{Total_Qualified}, \min(k, x))$
6. If $\text{Avg_Period}_{i-1} \neq 0$ $\text{Avg_Period}_i = \text{Avg_Period}_{i-1} * (\text{Qualified_Pulse_Period} / \text{Avg_Period}_{i-1})^k$
7. If $\text{Avg_Period}_{i-1} = 0$ $\text{Avg_Period}_i = \text{Qualified_Pulse_Period}$
8. $\text{Rate} = (60 / \Delta t) / \text{Avg_Period}_i$
- 30 9. $\text{Rate_Age} = \text{Rate_Age} + k * (\text{Qualified_Period_Age} - \text{Rate_Age})$

where:

r_i is the pulse rate corresponding to *Qualified_Pulse_Period*, in BPM

the $i-1$ subscript denotes the previous qualified pulse.

Δt is the oximetry algorithm's sample interval in seconds

5 $60 / \Delta t$ is the number of samples per minute

x is a filter weight that targets a 7-second response time for typical adult pulse rates.

k is the final filter weight, based on both x and the differences between consecutive values of r_i . During the first few pulses, k is increased to at least $1 / Total_Qualified$ so that the initial qualified pulses will be weighted equally.

10 [0068] *Consecutive_Qualified* is the number of consecutive qualified pulses, and *Total_Qualified* is the total number of pulses qualified since the subsystem was reinitialized. Both *Consecutive_Qualified* and *Total_Qualified* are incremented each time a *Qualified_Pulse_Period* is used, before k is calculated. *Consecutive_Qualified* is set to zero when a pulse is rejected by the pulse-period source currently in use.

15 [0069] The update formula for *Avg_Period_i*, in step 6 above, is a geometric average of *Avg_Period_i* and *Qualified_Pulse_Period*. Geometric averaging helps to keep the subsystem responsive to large pulses-to-pulse period variations, and large, sustained changes in pulse rate.

20 [0070] Once *Rate* is initialized to a non-zero value, *Rate_Age* is incremented every sample, whether or not *Rate* is updated.

Context Diagram

25 [0071] Fig. 3 is a context diagram of the pulse rate calculation subsystem. The subsystem updates its *Avg_Period* and *Rate* outputs from *Qualified_Pulse_Periods*. It uses *Qualified_Pulse_Periods* from the *Alternative_Period_Source* only if it last received a *Notify_Pulse_Rejected* from the primary source. It updates its *Rate_Age* output based on *Qualified_Period_Age*. When *Rate* is updated, the subsystem sets its *Pulse_Rate_Updated* flag. The *Reinitialize* input tells the subsystem to reinitialize itself. *Increment_Rate_Age* notifies the subsystem to increment its *Rate_Age* every sample once *Rate* is initialized.